

Product Datasheet

Helicobacter pylori Urease subunit beta antibody (orb240675)

Description Helicobacter pylori Urease subunit beta antibody

Species/Host Rabbit

Reactivity H. pylori

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Recombinant Helicobacter pylori Urease subunit beta protein (1-569AA)

Target ureB

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **P69996**

Expiration Date 12 months from date of receipt.



Western blot analysis of 293T whole cell...