

## Product Datasheet

# HPV18 Major capsid protein L1 antibody (FITC) (orb240689)

|                        |  |
|------------------------|--|
| <b>Description</b>     | HPV18 Major capsid protein L1 antibody (FITC)                  |
| <b>Species/Host</b>    | Rabbit   |
| <b>Reactivity</b>      | Virus  |
| <b>Conjugation</b>     | FITC   |
| <b>Immunogen</b>       | Recombinant HPV18 Major capsid protein L1 protein (1-568AA)    |
| <b>Target</b>          | L1   |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol |
| <b>Form/Appearance</b> | Liquid   |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | IgG  |
| <b>Clonality</b>       | Polyclonal   |
| <b>Uniprot ID</b>      | <b>P06794</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.                                |