


Product Datasheet

40S ribosomal protein S14 antibody (FITC) (orb240692)

Description	40S ribosomal protein S14 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human 40S ribosomal protein S14 protein (2-146AA)
Target	RPS14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62263
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...