


Product Datasheet

HE4 antibody (orb240758)


Description	HE4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human WAP four-disulfide core domain protein 2 protein (31-124AA)
Target	WFDC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14508
Dilution Range	WB:1:500-1:2000, IHC-P:1:200-1:500, IF/ICC:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of hepG2 cell



mouse heart using HE4 antibody



Immunohistochemical staining of human ov...