

Product Datasheet

JC polyomavirus Minor capsid protein VP2 antibody (Biotin) (orb240771)

Description	JC polyomavirus Minor capsid protein VP2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA)
Target	Minor capsid protein VP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03095
Expiration Date	12 months from date of receipt.