

## Product Datasheet

### JC polyomavirus Minor capsid protein VP2 antibody (FITC) (orb240772)

|                        |                                                                        |
|------------------------|------------------------------------------------------------------------|
| <b>Description</b>     | JC polyomavirus Minor capsid protein VP2 antibody (FITC)               |
| <b>Species/Host</b>    | Rabbit                                                                 |
| <b>Reactivity</b>      | Virus                                                                  |
| <b>Conjugation</b>     | FITC                                                                   |
| <b>Immunogen</b>       | Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA) |
| <b>Target</b>          | Minor capsid protein VP2                                               |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol         |
| <b>Form/Appearance</b> | Liquid                                                                 |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.          |
| <b>Note</b>            | For research use only                                                  |
| <b>Isotype</b>         | IgG                                                                    |
| <b>Clonality</b>       | Polyclonal                                                             |
| <b>Uniprot ID</b>      | <b>P03095</b>                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                        |