

Product Datasheet

Lactoplanarum L-lactate dehydrogenase 1 antibody (Biotin) (orb240779)

Description	Lactoplanarum L-lactate dehydrogenase 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Lactobacillus plantarum L-lactate dehydrogenase 1 protein (1-320AA)
Target	ldh1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56512
Expiration Date	12 months from date of receipt.