

Product Datasheet

Lemna minor ATP synthase subunit beta, chloroplastic antibody (FITC) (orb240784)

Description	Lemna minor ATP synthase subunit beta, chloroplastic antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Lemna minor ATP synthase subunit beta, chloroplastic protein (1-497AA)
Target	atpB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A9L9A3
Expiration Date	12 months from date of receipt.