

Product Datasheet

Lilrb3 Antibody, Biotin conjugated (orb240787)

| | |
|----------------------------|---|
| Description | Lilrb3 Antibody, Biotin conjugated |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Mouse Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (664-841AA) |
| Target | Lilrb3 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P97484 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)