

Product Datasheet

Macaca fascicularis IgG R FcRn large subunit p51 antibody (Biotin) (orb240791)

Description	Macaca fascicularis IgG R FcRn large subunit p51 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Macaca fascicularis IgG receptor FcRn large subunit p51 protein (24-297AA)
Target	FCGRT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8SPV9
Expiration Date	12 months from date of receipt.