

Product Datasheet

Moloney murine leukemia virus Gag polyprotein antibody (Biotin) (orb240795)

Description	Moloney murine leukemia virus Gag polyprotein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Moloney murine leukemia virus Gag polyprotein protein (216-478AA)
Target	gag
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03332
Expiration Date	12 months from date of receipt.