

Product Datasheet

C-C motif chemokine 12 antibody (Biotin) (orb240871)

Description	C-C motif chemokine 12 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse C-C motif chemokine 12 protein (23-104AA)
Target	Ccl12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q62401
Expiration Date	12 months from date of receipt.