

Product Datasheet

Hyaluronan synthase 2 antibody (FITC) (orb240982)

| | |
|------------------------|--|
| Description | Hyaluronan synthase 2 antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | FITC |
| Immunogen | Recombinant Mouse Hyaluronan synthase 2 protein (67-374AA) |
| Target | Has2 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P70312 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |