

Product Datasheet

VEGF antibody (Biotin) (orb241157)

Description	VEGF antibody (Biotin)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Pig Vascular endothelial growth factor A protein (27-190AA)
Target	VEGFA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49151
Expiration Date	12 months from date of receipt.