

Product Datasheet

Rabies virus Nucleoprotein antibody (FITC) (orb241176)

Description	Rabies virus Nucleoprotein antibody (FITC)
Species/Host	Rabbit
Reactivity	Rabies virus
Conjugation	FITC
Immunogen	Recombinant Rabies virus Nucleoprotein protein (1-450AA)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15197
Expiration Date	12 months from date of receipt.