

Product Datasheet

Rabies virus Phosphoprotein antibody (Biotin) (orb241181)

Description	Rabies virus Phosphoprotein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rabies virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rabies virus Phosphoprotein protein (1-297AA)
Target	P
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P69480
Expiration Date	12 months from date of receipt.