

Product Datasheet

Toxoplasma gondii Inflammatory profilin antibody (FITC) (orb241400)

Description	Toxoplasma gondii Inflammatory profilin antibody (FITC)
Species/Host	Rabbit
Reactivity	Parasite
Conjugation	FITC
Immunogen	Recombinant Toxoplasma gondii Profilin protein (2- 163AA)
Target	PRF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q58NA1
Expiration Date	12 months from date of receipt.