


Product Datasheet

Helicobacter pylori Urease subunit alpha antibody (orb241451)

Description	Helicobacter pylori Urease subunit alpha antibody
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Helicobacter pylori Urease subunit alpha protein (1-238AA)
Target	ureA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14916
Expiration Date	12 months from date of receipt.



Western blot analysis of Helicobacter py...