


Product Datasheet

Actin-related protein 2/3 complex subunit 2 antibody (Biotin) (orb241476)

Description	Actin-related protein 2/3 complex subunit 2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 2 protein (1-250AA)
Target	ARPC2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15144
Expiration Date	12 months from date of receipt.



Western blot analysis of using Actin-rel...