



Product Datasheet

Annexin A1 antibody (orb241487)


Description	Annexin A1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant Human Annexin A1 protein (5-346AA)
Target	ANXA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04083
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...



Immunofluorescence analysis of HeLa cell...



Immunoprecipitation analysis of HeLa who...