

Product Datasheet

Keratin, type II cytoskeletal 4 antibody (orb241679)

Description Keratin, type II cytoskeletal 4 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC

Immunogen Recombinant Human Keratin, type II cytoskeletal 4 protein (1-534AA)

Target KRT4

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P19013**

Dilution Range IHC-P:1:20-200, IF:1:50-200


Expiration Date 12 months from date of receipt.



Immunofluorescent analysis of A549 cells...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...