



Product Datasheet

Forkhead box protein P1 antibody (orb241937)


Description	Forkhead box protein P1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC
Immunogen	Recombinant Human Forkhead box protein P1 protein (1-114AA)
Target	FOXP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H334
Dilution Range	WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (1.1*...



Western blot analysis of Jurkat whole c...



Immunofluorescence staining of HeLa cell...