

Product Datasheet

Oxysterols R LXR-alpha antibody (Biotin) (orb241995)

Description	Oxysterols R LXR-alpha antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Oxysterols receptor LXR-alpha protein (1-402AA)
Target	NR1H3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13133
Expiration Date	12 months from date of receipt.