

Product Datasheet

R tyrosine-protein kinase erbB-2 antibody (Biotin) (orb242026)

Description	R tyrosine-protein kinase erbB-2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Receptor tyrosine-protein kinase erbB-2 protein (153-598AA)
Target	ERBB2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.