



Product Datasheet

Corticosteroid-binding globulin antibody (orb242117)


Description	Corticosteroid-binding globulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Corticosteroid-binding globulin protein (23-405AA)
Target	SERPINA6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08185
Dilution Range	WB:1:1000-1:5000, IHC-P:1:500-1:1000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Rat liver tiss...



Western blot aanalysis of U87 whole cell...