

Product Datasheet

Cytochrome b-c1 complex subunit Rieske, mitochondrial antibody (Biotin) (orb242186)

Description	Cytochrome b-c1 complex subunit Rieske, mitochondrial antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cytochrome b-c1 complex subunit Rieske, mitochondrial protein (79-274AA)
Target	UQCRFS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P47985
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...