



Product Datasheet

Fatty acid-binding protein antibody (orb242206)

Description	Fatty acid-binding protein antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Fatty acid-binding protein, liver protein (1-127AA)
Target	FABP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07148
Dilution Range	WB:1:500-2000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of recombinant Fa...



Immunohistochemistry of paraffin-embedde...