



Product Datasheet

Gamma-aminobutyric acid R-associated protein-like 2 antibody (orb242213)


Description	Gamma-aminobutyric acid R-associated protein-like 2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor-associated protein-like 2 protein (1-117AA)
Target	GABARAPL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60520
Dilution Range	WB:1:500-1:5000, IHC-P:1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...



Immunohistochemical staining of human heart...



Immunohistochemical staining of human prostate...