


Product Datasheet

Gamma-aminobutyric acid R-associated protein-like 2 antibody (Biotin) (orb242215)

Description	Gamma-aminobutyric acid R-associated protein-like 2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor-associated protein-like 2 protein (1-117AA)
Target	GABARAPL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60520
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...