



Product Datasheet

B-cell R-associated protein 31 antibody (orb242366)


Description	B-cell R-associated protein 31 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human B-cell receptor-associated protein 31 protein (2-243AA)
Target	BCAP31
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51572
Dilution Range	WB:1:1000-1:5000, IHC-P:1:1000-1:2000, IF:1:200-1:500, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...



Immunoprecipitation analysis of Mouse li...



Immunofluorescence analysis of HeLa cell...