

Product Datasheet

cytomegalovirus Envelope glycoprotein L antibody (FITC) (orb242369)

Description	cytomegalovirus Envelope glycoprotein L antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Human cytomegalovirus Envelope glycoprotein L protein (31-278AA)
Target	gL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	F5HCH8
Expiration Date	12 months from date of receipt.