


Product Datasheet

anemia virus apoptin antibody (Biotin) (orb242470)

Description	anemia virus apoptin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Chicken anemia virus Apoptin protein (1-121AA)
Target	VP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99152
Expiration Date	12 months from date of receipt.



Western blot analysis of Apoptin fusion ...