

Product Datasheet

Ovalbumin antibody (HRP) (orb242476)

Description	Ovalbumin antibody (HRP)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Chicken Ovalbumin (Native Protein)
Target	SERPINB14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01012
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.