

Product Datasheet

Hepatitis E virus genotype 1 Protein ORF3 antibody (Biotin) (orb242566)

Description	Hepatitis E virus genotype 1 Protein ORF3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis E virus genotype 1 Protein ORF3 protein (1-114AA)
Target	ORF3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O90299
Expiration Date	12 months from date of receipt.