

## Product Datasheet

# Hepatitis E virus genotype 1 Protein ORF3 antibody (FITC) (orb242567)

<b>Description</b>	Hepatitis E virus genotype 1 Protein ORF3 antibody (FITC)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	FITC
<b>Immunogen</b>	Recombinant Hepatitis E virus genotype 1 Protein ORF3 protein (1-114AA)
<b>Target</b>	ORF3
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">O90299</a>
<b>Expiration Date</b>	12 months from date of receipt.