

Product Datasheet

HIV-1 Gag-Pol polyprotein antibody (orb242569)

Description	HIV-1 Gag-Pol polyprotein antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag-Pol polyprotein protein (1-440AA)
Target	gag-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04585
Expiration Date	12 months from date of receipt.