

Product Datasheet

HIV-1 Gag-Pol polyprotein antibody (FITC) (orb242571)

Description	HIV-1 Gag-Pol polyprotein antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Human immunodeficiency virus type 1 group M subtype B Gag-Pol polyprotein protein (1-440AA)
Target	gag-pol
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04585
Expiration Date	12 months from date of receipt.