



## Product Datasheet

# Alanine aminotransferase 1 antibody (orb242594)


|                            |   |
|----------------------------|---|
| <b>Description</b>         | Alanine aminotransferase 1 antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IF, IHC, IP, WB  |
| <b>Immunogen</b>           | Recombinant Human Alanine aminotransferase 1 protein (2-496AA)  |
| <b>Target</b>              | GPT   |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <b>P24298</b>   |
| <b>Dilution Range</b>      | WB:1:500-2000, IHC-P:1:20-500, IF:1:50-200, IP:1:200-2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Rat heart tiss...



Western blot aanalysis of HepG2 whole ce...