

Product Datasheet

Ankyrin repeat and SOCS box protein 9 antibody (orb242614)

| | |
|----------------------------|---|
| Description | Ankyrin repeat and SOCS box protein 9 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human Ankyrin repeat and SOCS box protein 9 protein (1-294AA) |
| Target | ASB9 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q96DX5 |
| Expiration Date | 12 months from date of receipt. |