

Product Datasheet

herpesvirus Envelope glycoprotein D antibody (HRP) (orb242873)

Description	herpesvirus Envelope glycoprotein D antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human herpesvirus 1 Envelope glycoprotein D protein (26-339AA)
Target	gD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57083
Expiration Date	12 months from date of receipt.