



Product Datasheet

Interleukin-36 R antagonist antibody (orb242935)

Description	Interleukin-36 R antagonist antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Interleukin-36 receptor antagonist protein (1-155AA)
Target	IL36RN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBH0
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat kidney tis...



Western blot analysis of Mouse stomach ...