


Product Datasheet

MOB kinase activator 1B antibody (orb243020)


Description	MOB kinase activator 1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	Recombinant Human MOB kinase activator 1B protein (2-216AA)
Target	MOB1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7L9L4
Dilution Range	IHC-P:1:500-1:2000, IF:1:200-1:500, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of jurkat whole ce...



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human ap...