



Product Datasheet

Multiple inositol polyphosphate phosphatase 1 antibody (orb243028)

Description	Multiple inositol polyphosphate phosphatase 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Multiple inositol polyphosphate phosphatase 1 protein (31-487AA)
Target	MINPP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UNW1
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...



Immunohistochemical staining of human pl...