

Product Datasheet

papillomavirus type 18 Protein E7 antibody (FITC) (orb243084)

Description	papillomavirus type 18 Protein E7 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Human papillomavirus type 18 Protein E7 protein (1-105AA)
Target	E7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06788
Expiration Date	12 months from date of receipt.