



Product Datasheet

Probable global transcription activator SNF2L2 antibody (orb243110)


Description	Probable global transcription activator SNF2L2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC
Immunogen	Recombinant Human Probable global transcription activator SNF2L2 protein (700-1216AA)
Target	SMARCA2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51531
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...



Immunoprecipitation analysis of HeLa cel...



Immunohistochemical staining of human pa...