



Product Datasheet

Protein argonaute-2 antibody (orb243130)


Description	Protein argonaute-2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Recombinant Human Protein argonaute-2 protein (517-818AA)
Target	AGO2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UKV8
Dilution Range	WB:1:1000-1:5000, IHC-P:1:100-1:1500
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(Lane 1), J...



Immunohistochemical staining of human br...



Immunohistochemical staining of human ga...