


Product Datasheet

effector Noc2 antibody (orb243175)

Description	effector Noc2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Rab effector Noc2 protein (1-315AA)
Target	RPH3AL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UNE2
Expiration Date	12 months from date of receipt.



Western blot analysis of rat liver tissue...