

Product Datasheet

Ras-related protein Rab-11A antibody (FITC) (orb243192)

| | |
|------------------------|---|
| Description | Ras-related protein Rab-11A antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | FITC |
| Immunogen | Recombinant Human Ras-related protein Rab-11A protein (2-213AA) |
| Target | RAB11A |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P62491 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |