

Product Datasheet

rhinovirus A serotype 89 Genome polyprotein antibody (Biotin) (orb243238)

Description	rhinovirus A serotype 89 Genome polyprotein antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human rhinovirus A serotype 89 Genome polyprotein protein (575-866AA)
Target	Genome polyprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07210
Expiration Date	12 months from date of receipt.