
Product Datasheet


UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase 110 kDa subunit antibody (orb243357)

Description

UDP-N-acetylglucosamine-peptide N-



subunit antibody



Immunohistochemical staining of human br...

Species/Host

Rabbit

Reactivity

Human

Conjugation


Unconjugated

Tested Applications

ELISA, IF, IHC

Immunogen

Recombinant Human UDP-N-acetylglucosamine--peptide N-acetylglucosaminyltransferase 110 kDa subunit protein (606-1022AA)



Immunohistochemical staining of human pl...

Target

OGT

Preservatives


 Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance

Liquid

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.



Immunofluorescence analysis of Hela cell...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[O15294](#)
Dilution Range

IHC-P:1:20-200, IF:1:50-200

Expiration Date

12 months from date of receipt.